

2023



2023



exoneo

Inspired by anatomy

+33 (0)9 52 13 53 35

users@exoneo.eu

@exoneo.eu



@exoneo_eu

exoneo.eu



«Like a human foot»

Technical Sheet


"At EXONEO we devote our time and energy to one mission: To revolutionize the response to disability by drawing inspiration from the functioning of the human body.

We have developed a prosthetic foot: Upya, innovative and accessible to everyone in the world.

The Upya prosthetic foot reproduces the natural movement of the human foot while providing the dynamism to resume an active life, providing unprecedented comfort and rare stability."



Non-contractual images



With Upya, you can :

- Adapt your walking pace.
- Maintain a symmetry, step length, distance and walking time equivalent to a non-amputee.
- Cross most environmental barriers (slopes, stairs, uneven ground, etc.).

Other advantages of Upya :

- Smooth passage of the step.
- Energy restitution in two stages (blade + spring).
- Ability to adapt to a wide range of daily activities.
- Possibility of taking part in moderate recreational activities such as golf, cycling and hiking in leisure time.

Guarantee :

- Prosthetic foot : 2 years.
- Cosmetic cover : 6 months



ONE FOOT



MAIN FEATURES



PATIENT



WEATHER RESISTANT



CONTENT




*delivered in kit form.

PERFORMANCE : K1, K2



CHOOSE YOUR CONFIGURATION



P3 Patient high activity
P4
P5
P3 Patient low activity
P4
P5

+33 (0)9 52 13 53 35
users@exoneo.eu
@exoneo.eu



«Like a human foot»

UPYA



MOBILITY : K1, K2



CODE no 2793074

Upya is refunded in France

Materials

- Blades : Fiberglass
- Injected parts : PA 6.6 Polymer reinforced with 30% fiber glass
- Metallic parts : Aluminium FORTAL 7075
Stainless Steel
- Cushion heel in Polyurethane
- Cosmetic cover in EVA

Specifications

- Amputation level: Transfemoral, transtibial & bilateral
- Universal sizes, fits 25-28 cm (+29 cm accessory)
- Max. user weight : 110 kg
- Water resistant
- Reduced pressure in the socket
- 20° toe flexion
- Foot warranty: 24 months
- Foot shell warranty: 6 months

Upya reference

01C C -P

Chose side
D : Right
G : Left


Chose color
1 : Light skin
2 : Dark skin

Chose configuration
3, 4 or 5 : see chart ✱





*CHOOSE YOUR CONFIGURATION

Patient low activity



Patient high activity





CARBON FIBER

«Like a human foot»

UPYA TC



MOBILITY : K1, K2



CODE no 2793074

Upya is refunded in France

Materials

- Blades : Core layers in Fiberglass
External layers in Carbon Fiber
- Injected parts : PA 6.6 Polymer
reinforced with 30% fiber glass
- Metallic parts : Aluminium FORTAL 7075
Stainless Steel
- Cushion heel in Polyurethane
- Cosmetic cover in EVA

Specifications

- Amputation level: Transfemoral, transtibial & bilateral
- Universal sizes, fits 25-28 cm (+29 cm accessory)
- Max. user weight : 110 kg
- Water resistant
- Reduced pressure in the socket
- 20° toe flexion
- Foot warranty: 24 months
- Foot shell warranty: 6 months

Upya reference

01C C -P -TC

Chose side
D : Right
G : Left


Chose color
1 : Light skin
2 : Dark skin

Chose configuration
3, 4 or 5 : see chart *



*CHOOSE YOUR CONFIGURATION

Patient low activity



Patient high activity





«Like a human foot»

CARBON FIBER

FULL BLACK

UPYA TCFB



MOBILITY : K1, K2

Materials

- Blades : Core layers in Fiberglass
External layers in Carbon Fiber
- Injected parts : PA 6.6 Polymer
reinforced with 30% fiber glass
+ black gloss finish
- Metallic parts : Aluminium FORTAL 7075
Stainless Steel
- Cushion heel in Polyurethane
- Cosmetic cover in EVA

Specifications

- Amputation level: Transfemoral, transtibial & bilateral
- Universal sizes, fits 25-28 cm (+29 cm accessory)
- Max. user weight : 110 kg
- Water resistant
- Reduced pressure in the socket
- 20° toe flexion
- Foot warranty: 24 months
- Foot shell warranty: 6 months

01C C P TCFB

Chose side
D : Right
G : Left

Chose color
1 : Light skin
2 : Dark skin

Chose configuration
3, 4 or 5 : see chart *



CODE no 2793074

Upya is refunded in France



*CHOOSE YOUR CONFIGURATION

Patient low activity



Patient high activity

